

Otto Bock®

QUALITY FOR LIFE



Avantgarde Teen 2/VR

Sit right and be a part of the action!

Be a part of the action! The Teen 2 has what it takes to do so, for its abduction frame makes it light and manoeuvrable. Removable and swing-away footrests ensure easy transfer into and out of the wheelchair and individual leg support. The Teen 2 is highly variable in other ways, too. The rear part of the frame grows along with the user in 2-cm increments. Thus the knee angle remains constant and no additional parts are required. Three frame lengths and four variants for the front frame with various footplate options make the Teen 2 a solid basis for individual seating management. For sitting well is the best way to really get in on the action.

Otto Bock HealthCare GmbH

Max-Näder-Str. 15 · 37115 Duderstadt/Germany · Phone +49 5527 848-1590 · Fax +49 5527 848-1676 · reha-export@ottobock.com · www.ottobock.com

Avantgarde Teen 2/VR



Avantgarde Teen 2 in Detail


Standard Equipment:

- High-strength aluminium frame with single crossbrace
- Aluminium footrest with optimized locking function and a functional design
- Detachable footrests can be selected
- Back with different possible height settings as standard equipment (adjusts in 2.5 cm increments with anatomically adaptable back straps)
- Aluminium side panels
- Height- and angle-adjustable footrest
- Precise setting of the rear wheel position thanks to optimized adapter
- Front frame offset can be selected

Technical Data

Seat width	22–36 cm
Seat depth	24–40 cm
Back height	17.5–40 cm
Anterior seat height	40.5–49.5 cm
Posterior seat height	35–48 cm
Permissible user weight	90 kg
Weight	10.7 kg
Seat inclination	15°

Some of the models shown feature special accessories at additional charge. Information regarding delivery, appearance, measurements and weight corresponds to data available at the time of printing and is subject to change without notice.



Otto Bock®

QUALITY FOR LIFE