


REF.	10.BI.1214	10.BI.1216	10.BI.1216	10.BI.1217
Descrição Description Descripción	BIOMBO TRIPLO curvo, articulação nylon, c/ rodas / Curved folding triple screen panel, nylon articulation, with wheels / Paravent Triple courbé, articulation en nylon, avec roues / Biombo Tres cuerpos curvo, articulación de nylon, con ruedas	BIOMBO TRIPLO curvo, articulação nylon, c/ tacos / Curved folding triple screen panel, nylon articulation, with end caps / Paravent Triple courbé, articulation en nylon, avec embouts / Biombo Tres cuerpos curvo, articulación de nylon, con conteras	BIOMBO TRIPLO curvo, articulação nylon, c/ rodas / Curved folding triple screen panel, nylon articulation, with wheels / Paravent Triple courbé, articulation en nylon, avec roues / Biombo Tres cuerpos curvo, articulación de nylon, con ruedas	BIOMBO TRIPLO curvo, articulação nylon, c/ tacos / Curved folding triple screen panel, nylon articulation, with end caps / Paravent Triple courbé, articulation en nylon, avec embouts / Biombo Tres cuerpos curvo, articulación de nylon, con conteras
Material Artigo / Material Article Materiel Article / Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing Acier avec finition époxy / Acero con acabado Epoxy		Aço Inoxidável / Stainless steel Acier inoxydable / Acero inoxidable	
Dimensões Artigo / Dimensions Article Article Dimensions Dimensiones Artículo	2080x460x1840mm			
Peso Líquido / Net Weight Poids Net / Peso Neto	13kg	12,9kg	12,8kg	12,7kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	13,8kg	13,7kg	13,6kg	13,5kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,1186m <sup>3</sup>			



REF.	10.BI.1204	10.BI.1206	10.BI.1206	10.BI.1207
Descrição Description Descripción	BIOMBO Duplo curvo, articulação nylon, c/ rodas / Curved folding double screen panel, nylon articulation, with wheels / Paravent double courbé, articulation en nylon, avec roues / Biombo dos cuerpos curvo, articulación de nylon, con ruedas	BIOMBO Duplo curvo, articulação nylon, c/ tacos / Curved folding double screen panel, nylon articulation, with end caps / Paravent double courbé, articulation en nylon, avec embouts / Biombo dos cuerpos curvo, articulación de nylon, con conteras	BIOMBO Duplo curvo, articulação nylon, c/ rodas / Curved folding double screen panel, nylon articulation, with wheels / Paravent double courbé, articulation en nylon, avec roues / Biombo dos cuerpos curvo, articulación de nylon, con ruedas	BIOMBO Duplo curvo, articulação nylon, c/ tacos / Curved folding double screen panel, nylon articulation, with end caps / Paravent double courbé, articulation en nylon, avec embouts / Biombo dos cuerpos curvo, articulación de nylon, con conteras
Material Artigo / Material Article Materiel Article / Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing Acier avec finition époxy / Acero con acabado Epoxy		Aço Inoxidável / Stainless steel Acier inoxydable / Acero inoxidable	
Dimensões Artigo / Dimensions Article Article Dimensions Dimensiones Artículo	1405x460x1840mm			
Peso Líquido / Net Weight Poids Net / Peso Neto	9,4kg	9,3kg	9,2kg	9,1kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	10,2kg	10,1kg	10,0kg	9,9kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,1186m <sup>3</sup>			