



REF.	10.AD.1465	10.AD.1466
Descrição / Description Description / Descripción	Suporte Recolha de Sangue / Arm blood Collection Support / Support bras de prélèvement sanguin / Soporte Brazo Recogida de Sangre	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition epoxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	Ø600x745-1085mm	
Peso máximo admissível / Admissible maximum weight / Poids maximum admissible / Peso máximo admisible	20kg	
Peso Líquido / Net Weight Poids Net / Peso Neto	3kg	3,4kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	3,4kg	3,8kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,205m <sup>3</sup>	



REF.	10.CD.1326	10.CD.1325
Descrição / Description Description / Descripción	Cadeira p/ Recolha Sangue / Blood Collection Chair Chaise de collecte de sang / Silla para recoger sangre	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition epoxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	760x750x600mm	
Peso máximo admissível / Admissible maximum weight / Poids maximum admissible / Peso máximo admisible	120kg	
Peso Líquido / Net Weight Poids Net / Peso Neto	8,2kg	9kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	9,2kg	10kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,394m <sup>3</sup>	



REF.	10.AD.1451	10.AD.1452
Descrição / Description Description / Descripción	Suporte Bacia Simples / Simple Basin Support Support de cuvette simple / Porta palanganas simples	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition epoxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	Ø600x1000mm	
Peso Líquido / Net Weight Poids Net / Peso Neto	2,9kg	2,7kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	3,2kg	3kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,180m <sup>3</sup>	



REF.	10.AD.1453	10.AD.1454
Descrição / Description Description / Descripción	Suporte Bacia e Tabuleiro / Basin and Tray Support Porte cuvette et plateau / Soporte Cuenca y Tablero	
Material Artigo Material Article Materiel Article Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition epoxy / Acero con acabado Epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	Ø600x1004mm	
Peso Líquido / Net Weight Poids Net / Peso Neto	6,1kg	5,9kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	6,4kg	6,2kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,202m <sup>3</sup>	