




REF.	10.AD.1450	10.AD.1455
Descrição / Description / Descripción	Suporte Bacia Duplo / Double Basin Support / Support de cuvettes double / Soporte Cuenca Doble	
Material Artigo / Material Article / Matériel Article / Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável / Stainless steel / Acier inoxydable / Acero inoxidable
Dimensões Artigo / Dimensions Article / Article Dimensions / Dimensiones Artículo	670x600x1000mm	
Peso Líquido / Net Weight / Poids Net / Peso Neto	6,1kg	5,9kg
Peso Bruto / Gross Weight / Poids Brut / Peso Bruto	6,4kg	6,2kg
Volume Bruto / Gross Volume / Volume Brut / Volumen Bruto	0,202m <sup>3</sup>	



REF.	10.AD.1460	10.AD.1460.I	10.AD.1461	10.AD.1461.I
Descrição / Description / Descripción	Suporte Soro pintado 2 ganchos / Painted Serum Support 2 hooks / Support Sérum peint 2 crochets / Soporte Suero pintado 2 ganchos		Suporte Soro inox 2 ganchos / Stainless steel Serum Support 2 hooks / Support Sérum en acier inoxydable 2 crochets / Soporte Suero en acero inoxidable 2 ganchos	
Material Artigo / Material Article / Matériel Article / Material artículo	Aço com acabamento Epoxy e 2 ganchos em plástico / Steel with Epoxy Finishing and 2 plastic hooks / Acier avec finition époxy et 2 crochets plastique / Acero con acabado Epoxy y 2 ganchos plástico	Aço com acabamento Epoxy e 2 ganchos em inox / Steel with Epoxy Finishing and 2 stainless steel hooks / Acier avec finition époxy et 2 crochets en acier inoxydable / Acero con acabado Epoxy y 2 ganchos en acero inoxidable	Aço inoxidável e 2 ganchos em plástico / Stainless steel and 2 plastic hooks / Acier inoxydable et 2 crochets plastique / Acero inoxidable y 2 ganchos plástico	Aço inoxidável e 2 ganchos em inox / Steel with Epoxy Finishing and 2 stainless steel hooks / Acier avec finition époxy et 2 crochets en acier inoxydable / Acero con acabado Epoxy y 2 ganchos en acero inoxidable
Dimensões Artigo / Dimensions Article / Article Dimensions / Dimensiones Artículo	Ø600x1000-1780mm			
Peso máximo admissível por gancho / Hook Admissible maximum weight / Poids maximum admissible pour crochet / Peso máximo admissible por gancho	2kg			
Peso Líquido / Net Weight / Poids Net / Peso Neto	3,1kg	3,2kg	3,5kg	3,6kg
Peso Bruto / Gross Weight / Poids Brut / Peso Bruto	3,4kg	3,5kg	3,8kg	3,9kg
Volume Bruto / Gross Volume / Volume Brut / Volumen Bruto	0,098m <sup>3</sup>			



REF.	10.AD.1415	10.AD.1417	10.AD.1416	10.AD.1418
Descrição / Description / Descripción	Candeeiro Haste Flexível com lâmpada LED / Examination Lamp with flexible shaft and LED lamp / Lampe d'examen avec tige flexible et LED lampe / Candelabro hasta Flexible con lâmpara LED		Candeeiro Haste Flexível sem lâmpada LED / Examination Lamp with flexible shaft without LED lamp / Lampe d'examen avec tige flexible sans LED lampe / Candelabro hasta Flexible sin lâmpara LED	
Material Artigo / Material Article / Matériel Article / Material artículo	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável / Stainless steel / Acier inoxydable / Acero inoxidable	Aço com acabamento Epoxy / Steel with Epoxy Finishing / Acier avec finition époxy / Acero con acabado Epoxy	Aço Inoxidável / Stainless steel / Acier inoxydable / Acero inoxidable
Dimensões Artigo / Dimensions Article / Article Dimensions / Dimensiones Artículo	Ø600x1000-1750mm			
Peso Líquido / Net Weight / Poids Net / Peso Neto	2,20kg		2,15kg	
Peso Bruto / Gross Weight / Poids Brut / Peso Bruto	2,50kg		2,45kg	
Volume Bruto / Gross Volume / Volume Brut / Volumen Bruto	0,200m <sup>3</sup>		0,200m <sup>3</sup>	