



REF.	10.CR.1570	10.CR.1571
Descrição / Description Description / Descripción	Carro Transporte Roupa Limpa Transportation trolley for clean linen Gueridon de distribution de linge propre Carro dstribucion de ropa limpia	
Material Artigo Material Article Materiél Article Material artículo	Aço com acabamento epoxy Steel with epoxy finishing Acier avec finition époxy Acero con acabado epoxy	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	940x490x1020mm	
Peso Líquido / Net Weight Poids Net / Peso Neto	27,4kg	17kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	28,4kg	18kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,583m ³	



REF.	10.CR.1571.A.3156
Descrição Description Description Descripción	Carro Transporte Roupa Limpa Transportation trolley for clean linen Gueridon de distribution de linge propre Carro dstribucion de ropa limpia
Material Artigo Material Article Materiél Article Material artículo	Aço Inoxidável Stainless steel Acier inoxydable Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	800x500x900mm
Peso Líquido / Net Weight Poids Net / Peso Neto	12,7kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	13,7kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,444m ³



REF.	10.MC.7112
Descrição Description Description Descripción	Mesa de Gessos, em inox Stainless steel plaster trolley Gueridons pour plâtres en acier inoxydable Carro Escayola en acero inoxidable
Material Artigo / Material Article Materiél Article / Material artículo	Aço Inoxidável / Stainless steel Acier inoxydable / Acero inoxidable
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	500x1000x860mm
Peso Líquido / Net Weight Poids Net / Peso Neto	14,5kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	15,5kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,571m ³



REF.	10.CR.1500
Descrição Description Description Descripción	Carro Bala Oxigénio Trolley for Oxygen Bottle Chariot porte-bouteille d'oxygène Carro bala de oxigeno
Material Artigo Material Article Materiél Article Material artículo	Aço com acabamento epoxy Steel with epoxy finishing Acier avec finition époxy Acero con acabado epoxy
Dimensões Artigo / Dimensions Article Article Dimensions / Dimensiones Artículo	410x415x1470mm
Peso Líquido / Net Weight Poids Net / Peso Neto	9,5kg
Peso Bruto / Gross Weight Poids Brut / Peso Bruto	10,2kg
Volume Bruto / Gross Volume Volume Brut / Volumen Bruto	0,265m ³