

Hercules Prime Order form



Large rehab car seat for children, juveniles & adults in vehicles

**up to 80 kg / 176.4 lbs
body height 130 to 180 cm /
51.1 to 70.9"**

**tested according to
ECE R 44/04**



**Order e-mail:
info.export@thomashilfen.de**

**Status: 01.01.2024
valid until further notice**

Address

Customer no.:

Company:

Street:

Postal Code:

City:

Commission:

Contact person:

Phone:

E-mail:

Date:

BASIC MODEL **item code**

Hercules Prime rehab car seat

Please select the side of the pre-assembled **belt guide**:
(Change of side possible at any time. – A second belt guide at the top is available as an accessory.)



right

OR



left

980000460000

Standard features:

Back with integrated height adjustment; seat with integrated seat depth adjustment; continuous back angle adjustment with quick release lever; breathable upholstery with flame-retardant cover fabric: head cushion, back cushion, pelvic cushion, seat cushion, crotch pad; 5-point belt including shoulder pads; belt guide for 3-point vehicle belt (can be used on the right or left); prepared mounting options for all accessories; removable covers, machine washable up to 40° C

Measurements: Seat depth: 40 - 46 cm / 15.7 - 18.1" // Seat width: 38 cm / 15" // Back height (seat): 63 - 76 cm / 24.8 - 29.9"

REQUIRED OPTION FOR STABILITY: Always use the positioning seat with Seatfix adapter (419) or stabilizing bow (421).

Please chose ONE option:

+ **Seatfix adapter, behind the seat** * (if vehicle is equipped with Isofix)

9800004190000

+ **Stabilizing bow** * (if the vehicle is NOT equipped with Isofix)

9800004210000

ACCESSORIES **item code**

When function ordered in combination with basic model:

Please always select whether the desired function should or should not already be mounted on the basic model.

Swivel function

9800004010000

Swivel function PRE-MOUNTED

9800004011000

Please select the side (in driving direction) of the pre-assembled **swivel function**:
(Change of side possible)



right

OR



left

Options for further stability:

+ **Stand** **

9800004220000

+ **Cushion for function base** ***

9800005130000

Recommendation: For a better support use the swivel function with stand (422) and/or cushion (513), if necessary.

Tilting function 10° (0° / 2,5° / 5° / 7,5° / 10°)

9800004020000

Tilting function 10° (0° / 2,5° / 5° / 7,5° / 10°) **PRE-MOUNTED**

9800004021000

Tilting function 12,5 - 20° (only possible in combination with the tilting function 10°.)

9800004770000

Options for further stability:

+ **Stand** **

9800004220000

+ **Cushion for function base** ***

9800005130000

Recommendation: For a better support use the tilting function with stand (422) and/or cushion (513), if necessary.

please find further accessories on page 2

Hercules Prime Order form

Customer no.:

Company

Contact person:

Status: 01.01.2024
valid until further notice

Order e-mail:
info.export@thomashilfen.de

ACCESSORIES (continued)	item code
<input type="checkbox"/> Headrest, height adjustable	9800004040000
<input type="checkbox"/> Shoulder guide	9800004270000
<input type="checkbox"/> Upper body guide, firm ****	9800004280000
<input type="checkbox"/> Upper body guide, swing-away ****	9800004061000
<input type="checkbox"/> Lateral trunk supports, firm ****	9800004290000
<input type="checkbox"/> Lateral trunk supports, swing-away ****	9800004051000
<input type="checkbox"/> Abduction guide	9800004121000
<input type="checkbox"/> Leg guide	9800004150000
<input type="checkbox"/> Footrest, foldable ** (height and angle adjustable)	9800004170001
<input type="checkbox"/> Chest clip	9800004230000
<input type="checkbox"/> Second belt guide, top	9800004250000
<input type="checkbox"/> Car seat protection mat incl. back	9800008141000

NOTE:

* Seatfix adapter and stabilizing bow cannot be used in combination

** Stand and footrest cannot be used in combination

*** Cushion recommended for vehicle seats with seat recess to compensate for the contact surface of the function plate and vehicle seat

**** Upper body guide firm, upper body guide swing-away, lateral trunk supports firm and lateral trunk supports swing-away cannot be used in combination



Swivel function



Tilting function 10°
(0° / 2,5° / 5° /
7,5° / 10°)




Extension Tilting
function 12,5 - 20°
(12,5° / 15° / 17,5° / 20°)



Cushion
for function base



Seatfix adapter,
behind the seat



Stabilizing bow



Headrest,
height adjustable



Shoulder guide



Upper body guide,
firm



Upper body guide,
swing-away



Lateral trunk
supports,
firm



Lateral trunk
supports,
swing-away



Abduction guide




Leg guide



Stand



Footrest, foldable
(height and angle
adjustable)



Chest clip



Second belt guide,
top



Car seat protection
mat incl. back